

Deciphering Epigenetic Landscapes: A Comparative Analysis of DNA Methylation in 2D versus 3D Cultured Colon Cancer Cells using Scale Bioscience's Single Cell Methylation Kit

Sanika Khare¹, Lin Lin¹, Ashley Woodfin¹, Hosu Sin¹, Jerushah Thomas¹, Dmitry Pokholok¹, Eric Pu^{1,2}, Maggie Nakamoto¹

¹Scale Bioscience, 3210 Merryfield Row, San Diego, California, 92121


²email: eric.pu@scale.bio

ABSTRACT

We utilized our Single Cell Methylation Kit to generate DNA libraries from colon cancer cells, comparing the traditional monolayer (2D) setup with the more dynamic spheroid (3D) configuration. The increasing interest in 3D cell culture, especially spheroids, is driven by its potential to better approximate in vivo conditions compared to the more conventional 2D cultures. Our spheroids were cultivated in ultra-low adherence plates, while the 2D cells were maintained in TC treated plates.


Our workflow was highly efficient, yielding significant cytosine coverage and sequencing saturation. Using UMAP analysis, we observed distinct clustering patterns for the 2D and 3D cultured cells. There was some overlap, which we attribute to the cells having been in the 3D environment for only 5 days, suggesting they retained characteristics from their 2D origins. Additionally, by examining Differentially Methylated Regions (DMR), we could gain insights into how cellular behavior and expression might differ based on culture conditions. Such an approach not only has implications for potential gene discovery but also provides deeper insights into the epigenetic nuances of cells based on their growth environment. In conclusion, merging our Single Cell Methylation Kit with 3D spheroid culture techniques has provided valuable insights. This approach might prove pivotal in enhancing the precision and understanding of DNA methylation sites and cellular heterogeneity.

Figure 1: Combinatorial indexing technology can be used to increase single cell sequencing throughput, while decreasing cost



Combinatorial indexing uses repeated rounds of distributing cells across plates and subsequent well-based barcoding. Although multiple cells receive the same barcode in each round, repeated rounds creates a unique combinatorial barcode for the molecules in each individual cell.

Figure 2: Multiplexing of 2D and 3D cultured HCT-116 cells in scMET-seq workflow




A. HCT-116 cells were cultured in regular cell culture dish for 2D culture vs. low attachment culture dish for 3D culture. Spheroids were cultured for 5 days (they form on day 1). Nuclei from 2D and 3D culture were then purified, fixed, and treated before evenly distributed onto TSM plates for 1st level indexing and multiplexing. **B.** Cell recovery rate from the initial starting materials. **C.** Passing cell counts and cell rate after sequencing.

CONCLUSIONS

- The ScaleBio™ scMET workflow generates clean data with low background
- The ScaleBio™ Single Cell Methylation Sequencing Kits provides high sensitivity and low background single cell methylation information. Data from 2D and 3D cultured HCT-116 cells shows that the kit can separate 2D and 3D cultured cells solely based on methylation information.

Figure 3: ScaleBio's scMET-seq kit produces high quality data with high mapping rate and low CH methylation percentage



A. Knee plot shows uniform signal, confident cell calling, and low background across 2D and 3D cultured HCT-116 cells. **B.** Metrics show good sequencing and mapping quality. **C.** The median CG methylation percentage is in line with human cells (~75%). **D.** The low median CH methylation percentage (< 0.5%) indicating complete conversion of non-methylated Cytosines in CH context.

Figure 4: UMAP clustering separated 2D and 3D cells based on single cell methylation sequencing results

