

Background

Problem

Although plate-based fixation-compatible workflows have increased accessibility of single-cell profiling these workflows can still be challenging to use with complex samples, such as tissue, due to cell loss during preparation of single cell suspensions, variable recovery throughout the workflow, and arduous sample preparation leading to prohibitively long workflows.

Solution

Here we increase usability and flexibility of the ScaleBio™ Single Cell RNA Sequencing Kit by adding cell hashing to the workflow. Addition of these hashes enables pooling of cells from dissociated tissue samples before any centrifugation steps, drastically increasing cell recovery in upstream steps and preserving more of the sample for downstream recovery. Additionally, these hashes can increase sample throughput beyond the current 96 samples enabling larger screens.

Workflow

Cells or nuclei were hashed and taken through the ScaleBio workflow. Signal and background of the hashes were examined in different workflows (pooling of the samples before and after washing) as well as in different sample types (cell lines, PBMCs, and nuclei isolated from dissociated mouse liver and kidney).

Figure 3: Barnyard experiments show low background and good ability to identify cells with the correct hash regardless of whether hashed samples are mixed before or after washing.




Figure 1: Combinatorial indexing technology uses a plate-based workflow to barcode cells, easily increasing throughput without the need for any instrumentation.




Figure 4: Multiplexing with individual and pooled hashes recovers expected proportions of PBMCs and does not impact RNA assay performance.




Figure 2: Cell hashing can be used to label and pool samples upstream of the ScaleBio scRNA workflow.





Figure 5: Cell hashing of tissue-derived nuclei yields quality data with good sample identification.



Conclusions

- ScalePlex hash oligos are robustly captured in the ScaleBio scRNA workflow, enabling higher sample throughput, introducing new stopping points, and simplifying workflows.
- Hashing shows high signal and low background when tested on barnyard samples, PBMCs, and tissue-derived nuclei despite pooling of hashed samples before any wash steps.
- Hash oligos can be used to efficiently trace cells back to their original sample in barnyard cells, PBMCs, and tissue-derived nuclei.