

Introduction


Single cell-omics has transformed our understanding of cellular heterogeneity and cell function. Technology and informatics advances have accelerated single cell analysis but limitations in sample and cell throughput and high library cost per cell remain, creating a barrier to conduct larger, more complex studies. Additionally, extensibility to new assay modalities like methylation has proven to be challenging.

Single cell combinatorial indexing MET-seq (sciMET) has previously been proposed to scale throughput and reduce cost of single-cell methylation profiling (e.g., Nichols et al. 2022). Nuclei pass through a unique combination of barcoded wells using pooling and splitting, resulting in a unique combination of barcodes being attached to DNA for each cell. The advantages of combinatorial indexing are: 1) cell fixation for sample storage and transport, 2) sample and cell scalability, and 3) sub-linear cost scaling resulting in lower cost per cell.


In this study, we use ScaleBio™ Single Cell Methylation Kit to investigate DNA methylation patterns during development and oncogenesis, focusing on methylation changes in prenatal and adult brain tissue, and IDH-mutant glioma tumor tissue. Our analysis of single-cell methylomes showcases high single cell recovery and robust cytosine coverage. We identified cell-type specific clusters during the developmental states using Differentially Methylated Regions (DMR) and hypo- and hyper methylated regions. Comparisons with known bulk methylation profiles revealed unique single-cell methylation profiles, unveiling cellular heterogeneity and trajectories obscured by (pseudo) bulk analysis.

In conclusion, the ScaleBio™ Single Cell Methylation Kit is a first of its kind commercial solution providing a comprehensive genome-wide view of single cell methylomes with scalable throughput, high sensitivity and specificity offering unprecedented insights into cellular heterogeneity and trajectories during developmental and disease-driven methylation.


scMET Analysis of Human Adult and Prenatal Brain Tissue Reveals Diverse Cell Methylomes



ScaleBio scMET kit: A Modular & Streamlined Workflow Generates Comprehensive & High Quality Data



scMET Analysis of Human IDH^{mut} Glioma Tumor Produces Distinct Cell Type Clusters Based on Methylomes



Conclusions

- The ScaleBio™ Single Cell Methylation Kit (scMET kit) provides the first commercial solution to analyze single cell methylation at scale
- The workflow generates clean data with highly sensitive chemistry and low background. Starting with nuclei isolated from cells or tissues the scMET kit can process up to 18,400 nuclei per experiment
- Demonstrated ability to analyze a wide range of inputs from PBMCs to complex samples like human brain and tumor tissues with robust cytosine coverage and expected CH and CG detection
- In conclusion, the ScaleBio™ scMET kit offers a robust and highly scalable analysis solution for genome-wide single cell methylation profiling